

Anti-Sirt4 Polyclonal Antibody For If/Icc, Elisa - Q9Y6E7

Item Number: CM0014188

Introduction

SIRT4 Rabbit Polyclonal Antibody (Q9Y6E7) validated for IF/ICC and ELISA. Targets human SIRT4, a mitochondrial sirtuin involved in metabolism, tumor suppression, and DNA damage response. Ideal for oncology and metabolic disease research.

[Learn More](#)

Parameter	Value
Product Name	SIRT4 Rabbit pAb
Remarks/Alias	SIR2L4; SIRT4
Species	Human
Gene ID (Human)	23409
Gene ID	23409
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human SIRT4 (NP_036372.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9Y6E7
Protein Weight	35kDa
Shipping	Ice bag