

Anti-Sirt5 Rabbit Polyclonal Antibody (Ko Validated) For Wb, Elisa - Q9Nxa8

Item Number: CM0029541



Introduction

KO-validated SIRT5 rabbit polyclonal antibody for WB and ELISA, reactive with human, mouse, and rat. Explore mitochondrial deacylase functions in metabolism and aging research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SIRT5 Rabbit pAb
Remarks/Alias	SIR2L5; T5
Species	Human
Gene ID (Human)	23408
Gene ID	23408
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 37-310 of human SIRT5 (NP_036373.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NXA8
Protein Weight	34kDa
Shipping	Ice bag