

Anti-Sirt5 (Sir2L5) Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nxa8

Item Number: CM0022854



Introduction

SIRT5 Rabbit pAb, validated for WB and ELISA, recognizes human SIRT5 (SIR2L5), a mitochondrial NAD-dependent desuccinylase and demalonylase involved in metabolism. Ideal for metabolic research.

[Learn More](#)

Parameter	Value
Product Name	SIRT5 Rabbit pAb
Remarks/Alias	SIR2L5; SIRT5
Species	Human
Gene ID	23408
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 216-310 of human SIRT5 (NP_036373.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NXA8
Protein Weight	34kDa
Shipping	Ice bag