


Anti-Sirt6 Monoclonal Antibody For Wb, Elisa - Q8N6T7

Item Number: CM0010072



Introduction

Rabbit monoclonal antibody against human SIRT6, validated for Western blot and ELISA. Suitable for DNA repair, aging, metabolic, and cancer research. Recognizes human SIRT6.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] SIRT6 Rabbit mAb
Remarks/Alias	SIR2L6; hSIRT6; SIRT6
Species	Human
GeneID (Human)	51548
GeneID	51548
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 256-355 of human SIRT6 (Q8N6T7).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8N6T7
Protein Weight	36 kDa/39 kDa
Shipping	Ice bag