


Anti-Sirt6 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8N6T7

Item Number: CM0021135



Introduction

SIRT6 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for DNA damage repair, metabolism, and aging research. SWISS: Q8N6T7.

[Learn More](#)

Parameter	Value
Product Name	SIRT6 Rabbit pAb
Remarks/Alias	SIR2L6; hSIRT6; SIRT6
Species	Human
GeneID (Human)	51548
GeneID	51548
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-355 of human SIRT6 (NP_057623.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N6T7
Protein Weight	36 kDa/39 kDa
Shipping	Ice bag