


Anti-Sirt6 Rabbit Monoclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - Q8N6T7

Item Number: CM0011169



Introduction

High-quality rabbit monoclonal antibody targeting human SIRT6 (SIR2L6/hSIRT6). Validated for WB, IP, IF/ICC, IHC-P, and ELISA. Cross-reacts with mouse and rat. Suitable for DNA repair, metabolism, and aging research.

[Learn More](#)

Parameter	Value
Product Name	SIRT6 Rabbit mAb
Remarks/Alias	SIR2L6; hSIRT6
Species	Human
GeneID (Human)	51548
GeneID	51548
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-355 of human SIRT6 (NP_057623.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N6T7
Protein Weight	36 kDa/39 kDa
Shipping	Ice bag