

## Anti-Sirt7 Rabbit Monoclonal Antibody For Wb, Elisa - Q9Nrc8

Item Number: CM0006815



### Introduction

Rabbit monoclonal anti-SIRT7 antibody validated for WB and ELISA. Detects human SIRT7 with cross-reactivity to mouse. Targets the NAD-dependent deacetylase involved in chromatin regulation and transcription. Ideal for epigenetic research.

[Learn More](#)

| Parameter        | Value                                                                               |
|------------------|-------------------------------------------------------------------------------------|
| Product Name     | SIRT7 Rabbit mAb                                                                    |
| Remarks/Alias    | SIR2L7; SIRT7                                                                       |
| Species          | Human                                                                               |
| Gene ID (Human)  | 51547                                                                               |
| Gene ID          | 51547                                                                               |
| Immunogen        | Synthetic peptide corresponding to amino acids 301-400 of human SIRT7 (NP_057622.1) |
| Source           | Rabbit                                                                              |
| Category         | Monoclonal Antibodies                                                               |
| Application      | WB, ELISA                                                                           |
| Cross Reactivity | Mouse                                                                               |
| SWISS            | Q9NRC8                                                                              |
| Protein Weight   | 45 kDa                                                                              |
| Shipping         | Ice bag                                                                             |