

Anti-Siva1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O15304

Item Number: CM0027097

Introduction

SIVA1 rabbit polyclonal antibody validated for WB, IHC-P, ELISA in human, mouse. Targets apoptosis regulatory protein Siva. Ideal for apoptosis and T-cell signaling studies.

[Learn More](#)

Parameter	Value
Product Name	SIVA1 Rabbit pAb
Remarks/Alias	SIVA; CD27BP; Siva-1; Siva-2; SIVA1
Species	Human
GeneID (human)	10572
GeneID	10572
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human SIVA1 (NP_068355.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	O15304
Protein Weight	19kDa
Shipping	Ice bag