

Anti-Six1 Rabbit Polyclonal Antibody For Wb And Elisa - Q15475

Item Number: CM0016591



Introduction

Rabbit polyclonal antibody against human SIX1 (homeobox protein SIX1). Validated for WB and ELISA; cross-reacts with human and mouse. Suitable for studying embryonic development, myogenesis, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	SIX1 Rabbit pAb
Remarks/Alias	BOS3, TIP39, DFNA23, SIX1
Species	Human
GeneID (Human)	6495
GeneID	6495
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-270 of human SIX1 (NP_005973.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15475
Protein Weight	32kDa
Shipping	Ice bag