

Anti-Six2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Npc8

Item Number: CM0027909

Introduction

Rabbit polyclonal antibody against SIX2, validated for WB and ELISA with human and rat reactivity. SIX2 is a transcription factor crucial for kidney, craniofacial, and stomach development.

[Learn More](#)

Parameter	Value
Product Name	SIX2 Rabbit pAb
Remarks/Alias	SIX2
Species	Human
Gene ID (Human)	10736
Gene ID	10736
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 120-174 of human SIX2 (NP_058628.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9NPC8
Protein Weight	32kDa
Shipping	Ice bag