

Anti-Six3 Rabbit Polyclonal Antibody For Wb And Elisa - O95343

Item Number: CM0018969

Introduction

Rabbit polyclonal antibody against human SIX3, suitable for Western blot and ELISA. Cross-reacts with mouse and rat. Useful for studying forebrain development, eye formation, and holoprosencephaly.

[Learn More](#)

Parameter	Value
Product Name	SIX3 Rabbit pAb
Remarks/Alias	HPE2; SIX3
Species	Human
Gene ID (Human)	6496
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95343
Protein Weight	35kDa
Shipping	Ice bag