

Anti-Six5 Rabbit Polyclonal Antibody For Wb, Elisa - Q8N196

Item Number: CM0021784

Introduction

Rabbit polyclonal antibody against human SIX5 (BOR2/DMAHP). Validated for Western blot and ELISA with mouse and rat cross-reactivity. Targets transcription factor involved in organogenesis and retinal development. Research use only.

[Learn More](#)

Parameter	Value
Product Name	SIX5 Rabbit pAb
Remarks/Alias	BOR2, DMAHP, SIX5
Species	Human
Gene ID (Human)	147912
Gene ID	147912
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-170 of human SIX5 (NP_787071.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q8N196
Protein Weight	75kDa
Shipping	Ice bag