

Anti-Skap2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - 075563

Item Number: CM0015261



Introduction

Rabbit polyclonal antibody against human SKAP2, a Src kinase-associated phosphoprotein involved in B-cell adhesion and signaling. Validated for WB, IF/ICC, ELISA applications and cross-reacts with human, mouse, and rat.

[Learn More](#)

Parameter	Value
Product Name	SKAP2 Rabbit pAb
Remarks/Alias	PRAP, RA70, SAPS, SCAP2, SKAP55R, SKAP-HOM, SKAP2
Species	Human
Gene ID (Human)	8935
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-120 of human SKAP2 (NP_003921.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075563
Protein Weight	41kDa
Shipping	Ice bag