

Anti-Skiv2L2 Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - P42285

Item Number: CM0026836



Introduction

Rabbit polyclonal antibody targeting human SKIV2L2 (MTR4/DOB1). Validated for WB, IHC-P, IP, ELISA; cross-reacts with mouse and rat. Investigate exosome RNA helicase in nuclear RNA surveillance and processing.

[Learn More](#)

Parameter	Value
Product Name	SKIV2L2 Rabbit pAb
Remarks/Alias	Dob1; Mtr4; SKIV2L2; fSAP118; KIAA0052
Species	Human
Gene ID	23517
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 1-304 of human SKIV2L2 (NP_056175.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P42285
Protein Weight	118kDa
Shipping	Ice bag