

Anti-Skp1 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - P63208

Item Number: CM0008099



Introduction

SKP1 Rabbit mAb for Western blot, IF/ICC, ELISA. Designed for human, mouse, rat SKP1 (P63208). Detects a 19kDa protein component of SCF ubiquitin ligase, essential for cell cycle and signal transduction research.

[Learn More](#)

Parameter	Value
Product Name	SKP1 Rabbit mAb
Remarks/Alias	OCP2, p19A, EMC19, SKP1A, OCP-II, TCEB1L, SKP1
Species	Human
Gene ID (Human)	6500
Gene ID	6500
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 64-163 of human SKP1 (P63208)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63208
Protein Weight	19kDa
Shipping	Ice bag