

Anti-Skp2 Monoclonal Antibody For Wb, Ihc-P, Elisa - Q13309

Item Number: CM0011607



Introduction

Recombinant anti-SKP2 rabbit monoclonal antibody, validated for Western blot, immunohistochemistry, and ELISA. Recognizes human SKP2 (Q13309), an F-box protein regulating cell cycle progression via ubiquitin-mediated degradation of p27 and other substrates.

[Learn More](#)

Parameter	Value
Product Name	SKP2 Rabbit mAb
Remarks/Alias	p45; FBL1; FLB1; FBXL1; SKP2
Species	Human
Gene ID	6502
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human SKP2 (Q13309)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q13309
Protein Weight	48 kDa
Shipping	Ice bag