

Anti-Skp2 Polyclonal Antibody For Wb, If/Icc, Elisa - Q13309

Item Number: CM0017469



Introduction

Rabbit polyclonal antibody targeting human SKP2 (p45, FBL1). Applications: WB, IF/ICC, ELISA. Predicted molecular weight: 48 kDa. Cross-reacts with human. For studying cell cycle regulation, ubiquitin proteasome pathway, and antiviral mechanisms.

[Learn More](#)

Parameter	Value
Product Name	SKP2 Rabbit pAb
Remarks/Alias	p45; FBL1; FLB1; FBXL1; SKP2
Species	Human
Gene ID	6502
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human SKP2 (NP_005974.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q13309
Protein Weight	48kDa
Shipping	Ice bag