

Anti-Skp2 Rabbit Polyclonal Antibody For Western Blot, If/Icc, Elisa - Q13309

Item Number: CM0027712



Introduction

Validated rabbit polyclonal antibody targeting SKP2 for Western blot, IF/ICC, and ELISA detection in human, mouse, and rat samples. Key tool for cell cycle and ubiquitination research.

[Learn More](#)

Parameter	Value
Product Name	SKP2 Rabbit pAb
Remarks/Alias	p45; FBL1; FLB1; FBXL1; SKP2
Species	Human
GeneID (Human)	6502
GeneID	6502
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SKP2 (NP_005974.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13309
Protein Weight	48kDa
Shipping	Ice bag