

Anti-Sla2 Rabbit Polyclonal Antibody For Wb, Elisa - Q9H6Q3

Item Number: CM0018528



Introduction

SLA2 Rabbit Polyclonal Antibody for WB, ELISA. Detects human SLA2 (Src-like-adapter 2), a negative regulator of T-cell receptor signaling. Suitable for immunology and signaling research.

[Learn More](#)

Parameter	Value
Product Name	SLA2 Rabbit pAb
Remarks/Alias	MARS; SLAP2; SLAP-2; C20orf156; SLA2
Species	Human
GeneID (Human)	84174
GeneID	84174
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human SLA2 (NP_115590.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9H6Q3
Protein Weight	29kDa
Shipping	Ice bag