

Anti-Slamf6 Rabbit Polyclonal Antibody For Wb, Elisa - Q96Du3

Item Number: CM0016436



Introduction

Rabbit polyclonal antibody targeting human SLAMF6 (CD352/NTB-A). Validated for WB and ELISA. SLAMF6 is a self-ligand receptor modulating immune cell activation and differentiation. Ideal for immunology research.

[Learn More](#)

Parameter	Specification
Product Name	SLAMF6 Rabbit pAb
Remarks/Alias	KALI; NTBA; CD352; KAL1b; Ly108; NTB-A; SF2000; SLAMF6
Species	Human
Gene ID (Human)	114836
Gene ID	114836
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-215 of human SLAMF6 (NP_001171643.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q96DU3
Protein Weight	37kDa
Shipping	Ice bag