

Anti-Slamf7 Polyclonal Antibody For Wb And Elisa - Q9Nq25

Item Number: CM0027885

Introduction

Anti-SLAMF7 polyclonal antibody (rabbit) for WB and ELISA, cross-reacts with human, mouse, rat. Targets SLAM family member 7, a receptor modulating immune cell activation. Useful for immunology research. 37kDa, UniProt Q9NQ25.

[Learn More](#)

Parameter	Value
Product Name	SLAMF7 Rabbit pAb
Remarks/Alias	19A, CS1, CD319, CRACC, SLAMF7
Species	Human
GeneID (Human)	57823
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-226 of human SLAMF7 (NP_067004.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQ25
Protein Weight	37kDa
Shipping	Ice bag