

Anti-Slamf7 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nq25

Item Number: CM0020185

Introduction

SLAMF7 rabbit polyclonal antibody suitable for WB and ELISA. Detects human and mouse SLAMF7, a SLAM family receptor involved in NK cell activation and immune regulation.

[Learn More](#)

Parameter	Value
Product Name	SLAMF7 Rabbit pAb
Remarks/Alias	19A; CS1; CD319; CRACC; SLAMF7
Species	Human
Gene ID (Human)	57823
Gene ID	57823
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-226 of human SLAMF7 (NP_067004.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NQ25
Protein Weight	37kDa
Shipping	Ice bag