

Anti-Slc12A1 Rabbit Polyclonal Antibody For Wb, Elisa - Q13621

Item Number: CM0013212

Introduction

Rabbit polyclonal antibody targeting human SLC12A1/NKCC2, raised against a synthetic peptide. Suitable for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for kidney ion transport research.

[Learn More](#)

Parameter	Specification
Product Name	SLC12A1 Rabbit pAb
Remarks/Alias	BSC; BSC1; CCC2; BSC-1; NKCC2; SLC12A1
Species	Human
GeneID (Human)	6557
GeneID	6557
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human SLC12A1 (NP_000329.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q13621
Protein Weight	121kDa
Shipping	Ice bag