


Anti-Slc12A2 Polyclonal Antibody For Wb, Ihc-P, Elisa - P55011

Item Number: CM0014621



Introduction

Rabbit anti-SLC12A2 polyclonal antibody validated for WB, IHC-P, ELISA. Detects human, mouse, rat NKCC1 (131 kDa), a basolateral Na-K-2Cl cotransporter crucial for cell volume and ion homeostasis.

[Learn More](#)

Parameter	Value
Product Name	SLC12A2 Rabbit pAb
Remarks/Alias	BSC; BSC2; CCC1; BSC-2; KILQS; NKCC1; PPP1R141; SLC12A2
Species	Human
GeneID (Human)	6558
GeneID	6558
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 120-280 of human SLC12A2 (NP_001037.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55011
Protein Weight	131kDa
Shipping	Ice bag