

Anti-Slc12A7 Polyclonal Antibody For Wb, Elisa - Q9Y666

Item Number: CM0017081



Introduction

SLC12A7 rabbit polyclonal antibody validated for WB and ELISA, reactive with human, mouse, rat. Targets K-Cl cotransporter KCC4 involved in ion homeostasis and hearing. Ideal for IVD assay development and biomedical research.

[Learn More](#)

Parameter	Value
Product Name	SLC12A7 Rabbit pAb
Remarks/Alias	KCC4; SLC12A7
Species	Human
GeneID (Human)	10723
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human SLC12A7 (NP_006589.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y666
Protein Weight	119kDa
Shipping	Ice bag