


Anti-Slc14A1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13336

Item Number: CM0007406



Introduction

Rabbit monoclonal antibody against human SLC14A1 (urea transporter 1, JK blood group antigen), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse. Ideal for studies of urea transport in erythrocytes and renal collecting duct.

[Learn More](#)

Parameter	Value
Product Name	SLC14A1 Rabbit mAb
Remarks/Alias	JK; UT1; UTE; HUT11; Jk(a); Jk(b); RACH1; RACH2; UT-B1; HUT11A; HsT1341; SLC14A1
Species	Human
GeneID (Human)	6563
GeneID	6563
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 289-389 of human SLC14A1. (NP_056949.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13336
Protein Weight	43kDa
Shipping	Ice bag