


Anti-Slc14A1/Urea Transporter 1 Polyclonal Antibody For Wb, Ihc-P, If-P, Elisa - Q13336

Item Number: CM0021465



Introduction

Rabbit polyclonal antibody against human SLC14A1 (urea transporter 1, JK blood group antigen). Validated for WB, IHC-P, IF-P, and ELISA in human, mouse, and rat samples. Suitable for studying urea transport and renal physiology.

[Learn More](#)

Parameter	Value
Product Name	SLC14A1 Rabbit pAb
Remarks/Alias	JK; UT1; UTE; HUT11; Jk(a); Jk(b); RACH1; RACH2; UT-B1; HUT11A; HsT1341; SLC14A1
Species	Human
GeneID (Human)	6563
GeneID	6563
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 158-257 of human SLC14A1 (NP_001295208.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13336
Protein Weight	43kDa
Shipping	Ice bag