


Anti-Slc16A2 Rabbit Polyclonal Antibody (Mct8) For Wb, Elisa - P36021

Item Number: CM0027074



Introduction

Rabbit polyclonal antibody targeting human SLC16A2/MCT8, validated for WB and ELISA. Reacts with human, mouse, and rat tissues, detecting a 60kDa band. Suitable for thyroid hormone transport and MCT8-deficiency research.

[Learn More](#)

Parameter	Value
Product Name	SLC16A2 Rabbit pAb
Remarks/Alias	AHDS; MCT7; MCT8; XPCT; MCT 7; MCT 8; MRX22; DXS128; DXS128E; SLC16A2
Species	Human
GeneID (Human)	6567
GeneID	6567
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36021
Protein Weight	60kDa
Shipping	Ice bag