

Anti-Slc16A3/Mct4 Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - O15427

Item Number: CM0009878

Introduction

CamelBio offers Anti-SLC16A3/MCT4 Rabbit Monoclonal Antibody (O15427) validated for WB, IF/ICC, IHC-P, ELISA. This antibody targets the monocarboxylate transporter 4, ideal for lactate metabolism and cancer metabolism studies.

[Learn More](#)

Parameter	Value
Product Name	SLC16A3/MCT4 Rabbit mAb
Remarks/Alias	MCT3; MCT4; MCT 3; MCT 4; MCT-3; MCT-4
Species	Human
GeneID (Human)	9123
GeneID	9123
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 410-465 of human SLC16A3/MCT4 (NP_004198.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human
SWISS	O15427
Protein Weight	49kDa
Shipping	Ice bag