

Anti-Slc16A4 Rabbit Polyclonal Antibody For Wb And Elisa - O15374

Item Number: CM0013200

Introduction

Rabbit polyclonal antibody against human SLC16A4 (MCT4/5), a probable monocarboxylate transporter. Validated for WB and ELISA with reactivity in human and rat. Suitable for membrane protein research.

[Learn More](#)

Parameter	Details
Product Name	SLC16A4 Rabbit pAb
Remarks/Alias	MCT4; MCT5; SLC16A4
Species	Human
Gene ID	9122 (Human)
Immunogen	Synthetic peptide corresponding to amino acids 388-487 of human SLC16A4 (NP_004687.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O15374
Protein Weight	54 kDa
Shipping	Ice bag