

## Anti-Slc16A4 Rabbit Polyclonal Antibody For Wb, Elisa - O15374

Item Number: CM0030399

### Introduction

Rabbit polyclonal antibody targeting human SLC16A4 (MCT4/MCT5), suitable for Western blot and ELISA. Cross-reacts with mouse. Supports monocarboxylate transporter studies.

[Learn More](#)

| Parameter        | Value                 |
|------------------|-----------------------|
| Product Name     | SLC16A4 Rabbit pAb    |
| Remarks/Alias    | MCT4; MCT5; SLC16A4   |
| Species          | Human                 |
| GeneID (Human)   | 9122                  |
| GeneID           | 9122                  |
| Immunogen        | Synthetic peptide     |
| Source           | Rabbit                |
| Category         | Polyclonal Antibodies |
| Application      | WB, ELISA             |
| Cross Reactivity | Human, Mouse          |
| SWISS            | O15374                |
| Protein Weight   | 54kDa                 |
| Shipping         | Ice bag               |