

Anti-Slc16A7 (Mct2) Polyclonal Antibody For Wb, If/Icc, Elisa - O60669

Item Number: CM0015439



Introduction

High-quality rabbit polyclonal antibody against human SLC16A7/MCT2, validated for WB, IF/ICC, and ELISA. Detects monocarboxylate transporter 2 in human and rat samples. Ideal for metabolic and cancer research.

[Learn More](#)

Parameter	Specification
Product Name	SLC16A7 Rabbit pAb
Remarks/Alias	MCT2; SLC16A7
Species	Human
GeneID (Human)	9194
GeneID	9194
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	O60669
Protein Weight	52kDa
Shipping	Ice bag