

Anti-Slc16A7 (Mct2) Rabbit Polyclonal Antibody For Wb, Elisa - O60669

Item Number: CM0026065



Introduction

Rabbit polyclonal antibody targeting human SLC16A7 (MCT2), suitable for WB and ELISA. Cross-reacts with rat. Validated for endogenous monocarboxylate transporter detection in metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	SLC16A7 Rabbit pAb
Remarks/Alias	MCT2; SLC16A7
Species	Human
Gene ID (Human)	9194
Immunogen	Synthetic peptide corresponding to amino acids 379-478 of human SLC16A7 (NP_001257551.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O60669
Protein Weight	52kDa
Shipping	Ice bag