


Anti-Slc17A5 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nra2

Item Number: CM0017063



Introduction

Rabbit polyclonal antibody against human SLC17A5 (Sialin), validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for research on lysosomal transport and sialic acid metabolism.

[Learn More](#)

Parameter	Value
Product Name	SLC17A5 Rabbit pAb
Remarks/Alias	SD; AST; NSD; SLD; ISSD; VEAT; SIASD; SIALIN; SLC17A5
Species	Human
GeneID	26503
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SLC17A5 (NP_036566.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9NRA2
Protein Weight	55kDa
Shipping	Ice bag