

## Anti-Slc17A6 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - Q9P2U8

Item Number: CM0016475

### Introduction

Rabbit polyclonal antibody targeting SLC17A6/VGLUT2, suitable for Western blot, immunoprecipitation-fluorescence, and ELISA. Cross-reacts with human, mouse, and rat. For synaptic vesicle and neurotransmitter transport research.

[Learn More](#)

| Parameter        | Value                                                                                                      |
|------------------|------------------------------------------------------------------------------------------------------------|
| Product Name     | SLC17A6 Rabbit pAb                                                                                         |
| Remarks/Alias    | DNPI; VGLUT2; SLC17A6                                                                                      |
| Species          | Human                                                                                                      |
| GeneID (Human)   | 57084                                                                                                      |
| GeneID           | 57084                                                                                                      |
| Immunogen        | A synthetic peptide corresponding to a sequence within amino acids 500-582 of human SLC17A6 (NP_065079.1). |
| Source           | Rabbit                                                                                                     |
| Category         | Polyclonal Antibodies                                                                                      |
| Application      | WB, IF-P, ELISA                                                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                                          |
| SWISS            | Q9P2U8                                                                                                     |
| Protein Weight   | 64 kDa                                                                                                     |
| Shipping         | Ice bag                                                                                                    |