

Anti-Slc17A7 Polyclonal Antibody For Wb, If/Icc, If-P, Elisa - Q9P2U7

Item Number: CM0014620

Introduction

High-quality anti-SLC17A7 rabbit polyclonal antibody validated for WB, IF/ICC, IF-P, ELISA in human and mouse. Targets vesicular glutamate transporter 1 (VGLUT1) for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	SLC17A7 Rabbit pAb
Remarks/Alias	BNPI, VGLUT1, SLC17A7
Species	Human
Gene ID	57030
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9P2U7
Protein Weight	62kDa
Shipping	Ice bag