


Anti-Slc19A1 (Reduced Folate Transporter) Rabbit Polyclonal Antibody For Wb, Elisa - P41440

Item Number: CM0012921



Introduction

Rabbit polyclonal antibody targeting human SLC19A1 (reduced folate transporter), suitable for Western blot and ELISA. Ideal for folate transport, methotrexate uptake, and cyclic dinucleotide immune signaling research.

[Learn More](#)

Parameter	Value
Product Name	SLC19A1 Rabbit pAb
Remarks/Alias	RFC; CHMD; FOLT; IFC1; REFC; RFC1; hRFC; IFC-1; MEGAF; RFT-1; hSLC19A1; SLC19A1
Species	Human
Gene ID (Human)	6573
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 452-591 of human SLC19A1 (NP_919231.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P41440
Protein Weight	65kDa
Shipping	Ice bag