


Anti-Slc1A4 Polyclonal Antibody For Wb, Elisa - P43007

Item Number: CM0016223



Introduction

Rabbit polyclonal antibody targeting SLC1A4 (ASCT1, SATT) validated for Western blot and ELISA. Reacts with human, mouse, and rat. Ideal for amino acid transport and neurological disorder research.

[Learn More](#)

Parameter	Specification
Product Name	SLC1A4 Rabbit pAb
Remarks/Alias	SATT, ASCT1, SPATCCM, SLC1A4
Species	Human
GeneID (human)	6509
GeneID	6509
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 140-220 of human SLC1A4 (NP_003029.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43007
Protein Weight	56kDa
Shipping	Ice bag