


Anti-Slc1A4 Polyclonal Antibody For Wb, Ihc-P, Elisa - P43007

Item Number: CM0025337



Introduction

Rabbit anti-SLC1A4 polyclonal antibody, targeting human neutral amino acid transporter ASCT1. Validated in WB, IHC-P, ELISA, cross-reacts with mouse and rat. Ideal for studying D-serine homeostasis and amino acid transport mechanisms.

[Learn More](#)

Parameter	Value
Product Name	SLC1A4 Rabbit pAb
Remarks/Alias	SATT; ASCT1; SPATCCM; SLC1A4
Species	Human
Gene ID	6509
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 400-500 of human SLC1A4 (P43007)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43007
Protein Weight	56kDa
Shipping	Ice bag