


Anti-Slc1A5 Rabbit Polyclonal Antibody For Wb And Elisa - Q15758

Item Number: CM0015377



Introduction

High-quality rabbit polyclonal antibody targeting human SLC1A5 (ASCT2), suitable for Western blot and ELISA, with cross-reactivity to rat. Ideal for studying glutamine transport, cancer metabolism, and viral receptor functions.

[Learn More](#)

Parameter	Value
Product Name	SLC1A5 Rabbit pAb
Remarks/Alias	R16; AAAT; ATBO; M7V1; RDRC; ASCT2; M7VS1; SLC1A5
Species	Human
GeneID	6510 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human SLC1A5 (NP_005619.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q15758
Protein Weight	57kDa
Shipping	Ice bag