


Anti-Slc1A5 Rabbit Polyclonal Antibody For Wb, Elisa - Q15758

Item Number: CM0028724



Introduction

High-affinity rabbit polyclonal antibody against human SLC1A5 (ASCT2), validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for neutral amino acid transport and metabolic signaling research.

[Learn More](#)

Parameter	Value
Product Name	SLC1A5 Rabbit pAb
Remarks/Alias	R16; AAAT; ATBO; M7V1; RDRC; ASCT2; M7VS1; SLC1A5
Species	Human
Gene ID (Human)	6510
Gene ID	6510
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15758
Protein Weight	57kDa
Shipping	Ice bag