

Anti-Slc1A6 Rabbit Polyclonal Antibody For Wb And Elisa - P48664

Item Number: CM0028024

Introduction

High-quality anti-SLC1A6 rabbit polyclonal antibody validated for Western blot and ELISA. Detects human and mouse EAAT4, a glutamate transporter essential in cerebellar Purkinje cells. Ideal for neuroscience and transport studies.

[Learn More](#)

Parameter	Value
Product Name	SLC1A6 Rabbit pAb
Remarks/Alias	EAAT4; SLC1A6
Species	Human
GenelD (human)	6511
GenelD	6511
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 155-265 of human SLC1A6 (NP_001259016.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P48664
Protein Weight	62kDa
Shipping	Ice bag