

Anti-Slc1A7 Rabbit Polyclonal Antibody For Wb And Elisa - O00341

Item Number: CM0014749

Introduction

High-quality rabbit polyclonal antibody against human SLC1A7 (EAAT5), suitable for Western blot and ELISA. Shows mouse cross-reactivity. Ideal for glutamate transporter and retinal research applications.

[Learn More](#)

Parameter	Value
Product Name	SLC1A7 Rabbit pAb
Remarks/Alias	AAAT; EAAT5; SLC1A7
Species	Human
GeneID (Human)	6512
GeneID	6512
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 480-560 of human SLC1A7 (NP_006662.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O00341
Protein Weight	61kDa
Shipping	Ice bag