

Anti-Slc1A7 Rabbit Polyclonal Antibody For Wb, Elisa - O00341

Item Number: CM0024042

Introduction

Rabbit polyclonal antibody against SLC1A7 (EAAT5), validated for WB and ELISA applications, with cross-reactivity to mouse. Suitable for retina and neuroscience research on glutamate transport and synaptic inhibition.

[Learn More](#)

Parameter	Value
Product Name	SLC1A7 Rabbit pAb
Remarks/Alias	SLC1A7; AAAT; EAAT5; solute carrier family 1 member 7
Species	Human
Gene ID (Human)	6512
Gene ID	6512
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 480-560 of human SLC1A7 (NP_006662.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	O00341
Protein Weight	17kDa/60kDa
Shipping	Ice bag