

Anti-Slc20A1 Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Wum9

Item Number: CM0024883

Introduction

Validated SLC20A1 rabbit polyclonal antibody for WB, IF/ICC, ELISA. Detects human, mouse, rat. Sodium-dependent phosphate transporter 1, relevant in calcification and retroviral receptor research.

[Learn More](#)

Parameter	Value
Product Name	SLC20A1 Rabbit pAb
Remarks/Alias	PIT1; GLVR1; PiT-1; Glvr-1; SLC20A1
Species	Human
Gene ID (Human)	6574
Gene ID	6574
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 257-356 of human SLC20A1 (NP_005406.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WUM9
Protein Weight	74kDa
Shipping	Ice bag