

Anti-Slc20A2 Polyclonal Antibody For Wb And Elisa - Q08357

Item Number: CM0028165

Introduction

Rabbit polyclonal antibody targeting SLC20A2 (PiT-2), validated for WB and ELISA with cross-reactivity to human, mouse, and rat. Suitable for research on phosphate transport, bone mineralization, and retroviral entry.

[Learn More](#)

Parameter	Value
Product Name	SLC20A2 Rabbit pAb
Remarks/Alias	PIT2; RAM1; GLVR2; IBGC1; IBGC2; IBGC3; MLVAR; PIT-2; Ram-1; GLVR-2; SLC20A2
Species	Human
GeneID (Human)	6575
GeneID	6575
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 235-485 of human SLC20A2 (NP_001244110.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q08357
Protein Weight	70kDa
Shipping	Ice bag