


Anti-Slc22A11 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nsa0

Item Number: CM0028454



Introduction

Anti-SLC22A11 Rabbit Polyclonal Antibody (OAT4) validated for WB, ELISA. Detects human, mouse, rat SLC22A11, a placental/kidney organic anion transporter. Ideal for transport studies and research on renal reabsorption and placental steroid exchange.

[Learn More](#)

Parameter	Value
Product Name	SLC22A11 Rabbit pAb
Remarks/Alias	OAT4; hOAT4; SLC22A11
Species	Human
GeneID (Human)	55867
GeneID	55867
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 40-140 of human SLC22A11 (NP_060954.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NSA0
Protein Weight	60kDa
Shipping	Ice bag