


Anti-Slc22A2 (Oct2) Polyclonal Antibody For Wb And Elisa - O15244

Item Number: CM0025390



Introduction

Rabbit polyclonal antibody targeting human SLC22A2/OCT2 organic cation transporter. Validated for Western blot and ELISA; cross-reacts with rat. Ideal for studying drug clearance, neurotransmitter transport, and polyamine pathways. Predicted molecular weight 63kDa.

[Learn More](#)

Parameter	Value
Product Name	SLC22A2 Rabbit pAb
Remarks/Alias	OCT2; SLC22A2
Species	Human
Gene ID (Human)	6582
Gene ID	6582
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 44-140 of human SLC22A2 (NP_003049.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O15244
Protein Weight	63kDa
Shipping	Ice bag