


## Anti-Slc22A3/Oct3 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O75751

Item Number: CM0003694



### Introduction

Recombinant rabbit monoclonal anti-SLC22A3/OCT3 antibody validated for WB, IHC-P, and ELISA. Detects human and mouse solute carrier family 22 member 3 (O75751), a key organic cation transporter involved in neurotransmitter and drug uptake.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | SLC22A3/OCT3 Rabbit mAb                                                                                              |
| Remarks/Alias    | EMT; EMT3; OCT3; SLC22A3/OCT3                                                                                        |
| Species          | Human                                                                                                                |
| GeneID (Human)   | 6581                                                                                                                 |
| GeneID           | 6581                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 301-450 of human SLC22A3/OCT3 (O75751) |
| Source           | Rabbit                                                                                                               |
| Category         | Monoclonal Antibodies                                                                                                |
| Application      | WB, IHC-P, ELISA                                                                                                     |
| Cross Reactivity | Human, Mouse                                                                                                         |
| SWISS            | O75751                                                                                                               |
| Protein Weight   | 61kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |