

Anti-Slc22A5 (Octn2) Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O76082

Item Number: CM0021871



Introduction

Rabbit polyclonal antibody against SLC22A5 (OCTN2), a high-affinity carnitine transporter. Validated for WB, IF/ICC, ELISA in human, mouse, rat. Suitable for studying carnitine metabolism, intestinal homeostasis, and primary carnitine deficiency.

[Learn More](#)

Parameter	Value
Product Name	SLC22A5 Rabbit pAb
Remarks/Alias	CDSP; OCTN2; SLC22A5
Species	Human
Gene ID (Human)	6584
Gene ID	6584
Immunogen	Recombinant fusion protein containing amino acids 430-529 of human SLC22A5 (NP_003051.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O76082
Protein Weight	63kDa
Shipping	Ice bag